

The *easiest* way to floss with braces!



Insert spatula end of the flosser under your wire and press against teeth.



With spatula against your teeth, slide floss between your teeth.



In difficult to reach areas, maintain pressure against your teeth while flossing.



The bracket brush cleans in and around your brace brackets.

#1
PROFESSIONALLY RECOMMENDED!

Platypus[®]Co.com
MADE IN AMERICA

Platypus[®]
Orthodontic
Flosser

Toll Free: 888.543.6313

NEW

Platypus[®] Orthodontic Toothbrush



Features:

- **SOFT OUTER BRISTLES** are angled inward to reach in and around brackets and arch-wire.
- **MEDIUM INNER BRISTLES** are shorter to target metal components for a more effective cleaning.



**DESIGNED SPECIFICALLY
FOR BRACES!**