ELECTRONIC INSTRUCTIONS FOR USE

PLATYPUS® ALIGNER SEATING CHEW – ENGLISH / EN

COMPANY IDENTIFICATION:

 Glenroe, LLC
555 W. Brown Deer Road, Suite 200
Milwaukee, WI 53217 USA

24 HR Emergency Phone: +1 406 626 1736
Phone / Information: +1 406 626 1736
Website: www.platypusco.com/product-downloads

PRODUCT INFORMATION

Name: Platypus® Aligner Seating Chew

 MT Promedt Consulting GmbH
Ernst-Heckel-Straße 7
66386 St. Ingbert, Germany

INTENDED PURPOSE / INDICATIONS FOR USE

Platypus Aligner Seating Chew is intended to be used to assist in the placement of an orthodontic aligner. The intended patient population of this device is any individuals undergoing orthodontic treatment with clear orthodontic aligners. The intended users of this device are orthodontists, dentists, and other clinically trained personnel involved in orthodontic treatment.

Platypus Aligner Seating Chew is indicated for use in patients undergoing orthodontic aligner treatment.

PERFORMANCE CHARACTERISTICS

The Platypus Aligner Seating Chew is a device that mechanically aids in the placement of an aligner.

CONTRAINDICATIONS & SIDE EFFECTS

- Do not swallow.
- Do not use on children without adult supervision.
- Do not use outside of manufacturer's instructions.
- Do not reuse.

WARNINGS

- Adult supervision is recommended for children.
- Tool should be used during aligner installation only as directed.
- Discard after use.

INSTRUCTIONS

GENERAL


1. Move aligner in place.
2. Position the aligner chew between teeth.
3. Begin biting on chew until aligner is firmly positioned on teeth.

STORAGE

No special storage conditions required.

DISPOSAL

Dispose of aligner seating chew and packaging material in accordance with local and governing laws.

STERILIZATION

This product is labeled as non-sterile. The devices are also designed, manufactured, and packaged as non-sterile.

SERIOUS INCIDENT REPORTING SECTION

If, during the use of this device or as a result of its use, a serious incident has occurred please report it to the manufacturer and/or the competent authority of the Member State in which the user and/or patient is established.

PAPER COPY REQUEST

A paper copy of this IFU can be provided upon request. To request, please reach out to the Glenroe team at regulatory@glenroe.com.

REGULATORY INFORMATION

	MANUFACTURER
	REF NUMBER
	USE-BY SYMBOL WITH EXPIRATION DATE
	NOT STERILIZED
 www.platypusco.com/product-downloads	ELECTRONIC INSTRUCTIONS FOR USE
	NOT MADE WITH NATURAL RUBBER LATEX
	MEDICAL DEVICE
	CE MARK
	LOT NUMBER
	EU REP

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